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PARI LC SPRINT® SINUS

Model: PARI LC SPRINT SINUS (Type 023)

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
Model: PARI LC SPRINT SINUS (Type 023)

Nebuliser for PARI SINUS / PARI SINUS N / PARI SINUS2 nebuliser systems

Important: Read these instructions carefully before using the product for the first time. Follow all instructions and safety instructions!

Keep the instructions in a safe place.

Identification, validity, version

These instructions for use are valid for PARI LC SPRINT SINUS nebuliser (Type 023) for PARI SINUS, PARI SINUS N and PARI SINUS2 nebuliser systems in the following countries:

Countries outside the EU

Version of these instructions for use: Version F – 2020-02, Approved version dated 2019-11-19

Information as of: 2019-11

The current version of the instructions for use can be downloaded from the internet as a PDF file:

www.pari.com (on the respective product page)

Formats available for visually impaired patients

The instructions for use available in PDF format on the internet can be enlarged for printing.

CE conformity

This product satisfies the requirements of 93/42/EEC (Medical devices).

Brands

The following brand names are registered trademarks of PARI GmbH Spezialisten für effektive Inhalation in Germany and/or other countries.

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No part of this documentation may be reproduced or processed, duplicated, translated or disseminated by electronic means without the prior, written consent of PARI GmbH.

All rights reserved. Technical and design changes and printing errors reserved. Principle shown in illustrations.

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Contact

For all product information and in the event of defects or questions about usage, please contact our Service Center:

Tel.: +49 (0)8151-279 220 (international)

+49 (0)8151-279 279 (German-speaking)

E-Mail: info@pari.de

TABLE OF CONTENTS

IMPORTANT INFORMATION	7
Intended purpose	7
Indication	8
Contraindications.....	8
Safety instructions	8
PRODUCT DESCRIPTION	11
Components	11
Product combinations	12
Description of function	12
Material information	12
Operating life	13
APPLICATION	13
Preparing for treatment	13
Performing treatment.....	16
Ending the treatment	20
REPROCESSING	21
Preparation	21
Cleaning	21
Disinfecting	22
Care of the PARI SINUS tubing system	23
Inspecting	23
Drying	24
Storage.....	24
TROUBLESHOOTING	24
TECHNICAL DATA	25
General nebuliser data	25
Aerosol data	25
MISCELLANEOUS	26
Disposal.....	26
Labelling	26

APPENDIX: Reprocessing in professional environments (for only one patient)	27
Nebuliser	27
PARI SINUS tubing system	29

1 IMPORTANT INFORMATION

1.1 Intended purpose

The PARI LC SPRINT SINUS nebuliser generates aerosols¹ for respiratory therapy.

Together with the PARI SINUS, PARI SINUS N or PARI SINUS2 compressor, the nebuliser forms a nebuliser system for treatment of the upper airways.

The nebuliser is designed to treat patients from age 6 up to adults.

Only solutions and suspensions that are approved for use in nebuliser therapy may be used.

The nebuliser must only be used with a PARI SINUS, PARI SINUS N or a PARI SINUS2 compressor.

For reasons of hygiene, this PARI product must only be used in a home environment by a single patient.

The nebuliser system must only be used by individuals who understand the contents of the instructions for use and are able to operate the nebuliser system safely. Individuals in the following groups must be supervised by a person who is responsible for their safety:

- Children
- Individuals with limited capabilities (e.g., physical, mental, sensory)

If the patient is not able to operate the nebuliser system safely himself, the therapy must be carried out by the responsible person.

An application takes between 5 and 10 minutes (depending on the quantity of fluid), but in no case more than 20 minutes.

The frequency and duration of use will be determined by the physician or therapist according to the individual needs of the patient.

1) Aerosol: Small particles of solid, liquid or mixed composition (fine "mist") suspended in gases or air.

1.2 Indication

Diseases of the nose and paranasal sinuses.

1.3 Contraindications

This product is only designed for patients who are able to breathe by themselves and are conscious.

Following surgery or where inflammation of the middle ear is present, a risk/benefit assessment for vibrating aerosols must be carried out by the doctor before the start of treatment.

1.4 Safety instructions

The present instructions for use contain important information, safety instructions and precautionary measures. The user must follow these in order to guarantee safe operation of this PARI product.

This PARI product must only be used as described in these instructions for use.

The instructions for use of the compressor and the information for use of the inhalation solution used must also be followed.

General

Nebuliser aerosol characteristics differ from the information provided by the manufacturer if non-approved solutions or suspensions are used for nebulisation.

This product is not suitable for use in an anaesthetic breathing system or a ventilator breathing system.

If your health condition is not improved, or even worsens as a result of the treatment, seek professional medical advice.²

2) Professional medical staff: Doctors, pharmacists and physiotherapists.

Hazard due to small parts which can be swallowed

The product contains small parts. Small parts can block the airways and lead to a choking hazard. Keep all components of the product out of the reach of babies and infants at all times.

Hygiene

Observe the following hygiene instructions:

- Do not use product components unless they have been thoroughly cleaned and dried. Contamination and residual moisture encourage the growth of bacteria, which increases the risk of infection.
- Before every use and cleaning, wash your hands thoroughly.
- It is essential to clean and dry all product components before the first application as well.
- Always use drinking water for cleaning and disinfecting.
- Dry all product components completely after each cleaning and disinfection.
- Do not keep the product components in a damp environment or together with damp objects.

Treatment of children and anyone who requires assistance

For individuals who are not able to perform the therapy session without assistance or cannot appreciate the hazards, the risk of injury is greater (e.g., strangulation with the power cord or a tube). Such individuals include for example babies, children and people with limited capabilities. For these individuals, a person responsible for safety must supervise or carry out the application.

Identifying and classifying warning instructions

Safety-critical warnings are categorised according to the following hazard levels in these instructions for use:

DANGER

DANGER indicates a hazardous situation which will lead to severe injuries or death if it is not avoided.

WARNING

WARNING indicates a hazardous situation which can lead to severe injuries or death if it is not avoided.

CAUTION

CAUTION indicates a hazardous situation which can lead to mild or moderate injuries if it is not avoided.


NOTE

NOTE indicates a situation which can lead to material damage if it is not avoided.

2 PRODUCT DESCRIPTION

2.1 Components

The following components are included in the package:

(1)	Nebuliser
(1a)	Nebuliser upper part
(1b)	Nozzle insert (orange)
(1c)	Nebuliser lower part
(1d)	Adapter for vibration
(1e)	Nasal plug ³
(1f)	Nasal joining piece
(1g)	Bend
(1h)	PARI SINUS tubing system ³
(2)	Air filter for compressor types 128/130/152 ³
(3)	Air filter for compressor type 028 ³

³) Not included with all product variants.


2.2 Product combinations

The PARI LC SPRINT SINUS nebuliser can only be operated with the PARI SINUS, PARI SINUS N and the PARI SINUS2 compressors.

2.3 Description of function

The nebuliser is part of a PARI SINUS nebuliser system.

When compressed air is supplied, the nebuliser generates an aerosol from the liquid, for example the medication with which it is filled. This aerosol is delivered to the nose and paranasal sinuses through the nasal insert. The pulsation supports this process.

2.4 Material information

The individual product components are made from the following materials:

Product component	Material
Nebuliser upper part	Polypropylene, thermoplastic elastomer
Nozzle insert	Polypropylene
Nebuliser lower part	Polypropylene, thermoplastic elastomer
PARI SINUS tubing system	Polyvinyl chloride
Tubing endpiece	Thermoplastic elastomer
Nasal joining piece	Polypropylene, thermoplastic elastomer
Adapter for vibration	Thermoplastic elastomer
Nasal plug	Silicone
Bend	Polypropylene

2.5 Operating life

The individual product components have the following expected operating lives:

Product component	Operating life
Nebuliser (all components except the PARI SINUS tubing system)	300 disinfections, max. 1 year
PARI SINUS tubing system	Max. 1 year

3 APPLICATION

People who assist others in carrying out the therapy must ensure that all of the steps described below are carried out correctly.

3.1 Preparing for treatment

Assembling the nebuliser


CAUTION

Risk of impaired treatment


Damaged components and/or an incorrectly assembled nebuliser may impair functioning of the nebuliser and thus treatment as well.

- Check all nebuliser components and the accessories before each use.
- Replace any broken, misshapen or seriously discoloured parts.
- Follow the assembly instructions in these instructions for use.


- Press the nozzle insert lightly onto the nozzle in the nebuliser lower part. The arrow on the nozzle insert must point upwards.
- Place the nebuliser upper part on the nebuliser lower part and turn it clockwise to lock the two parts together.

- Connect the nasal joining piece to the bend and then connect both to the nebuliser.

- Insert one tubing connector on the PARI SINUS tubing system in the bottom of the nebuliser.

Filling the nebuliser

NOTE


Nebuliser lid might break off

If the lid is twisted in the wrong direction, it may break off. The nebuliser will then be unusable and irreparable.

- Never move the lid except in the direction allowed by the hinge.

- Insert the nebuliser in the holder on the compressor intended for this purpose.
- Open the nebuliser lid by pressing your thumb against the underside of the lid.
- Pour the required quantity of inhalation solution into the top of the nebuliser.

Note the minimum and maximum fill volumes [see: General nebuliser data, page 25]. If the nebuliser contains too little or too much liquid, the nebulisation and consequently the therapy will be less effective.

- Close the nebuliser lid.
Make sure that the lid snaps into place.

If several inhalation solutions are to be used one after the other:


- Rinse the nebuliser out with drinking water between the individual applications.
- Shake excess water out of the nebuliser.
- Fill the nebuliser with the next inhalation solution as described.

Connecting the tubes

The PARI SINUS tubing system consists of two identical tubes. One tube is needed for the compressed air supply, the other helps to generate the vibration.


To connect the tubing system up fully, proceed as follows:

- Insert the corresponding connector of the unconnected tube in the vibration adapter.
- Insert the adapter in the nebuliser upper section.

- Insert the free tubing connectors in the corresponding connectors on the compressor.

Info: The tubing connectors can be inserted more easily if they are twisted slightly when inserting.

3.2 Performing treatment

All of the safety instructions in these instructions for use must have been read and understood before any treatment is carried out.

i Always hold the nebuliser upright during treatment.

General instructions for use

To ensure that the treatment has the desired effect, please follow the instructions below:

- In order for the aerosol to reach all regions of the nose and paranasal sinuses, the nose must be clear. If your nose is blocked, rinse it beforehand using a nasal douche for example. If necessary, use a product designed to reduce swelling.
- The aerosol that is generated during nebulisation should only be delivered to the nose and paranasal sinuses. Therefore, the soft palate⁴ – the natural partition between the nasal cavity and the mouth and throat – must be closed.

Tips for closing the soft palate

To close your soft palate, say the three letters “ing” as one word for as long as possible. After a short breathing pause, repeat this as often as necessary until the treatment is complete.

i *When you make certain sound combinations such as “ing” and when you swallow, the “soft palate” closes and separates the nasal cavity from the mouth and throat automatically.*

Alternatively, you can breathe in and out slowly and deliberately through your mouth, as if attempting to blow a candle out (i.e. do not breathe through your nose). However, this method of closing the soft palate is not easy to perform for everybody, and some practice may be required.

While the soft palate is closed and you are holding your breath, the aerosol flows into the paranasal sinuses. Your nostrils will vibrate a little.

i *Caution: If the soft palate is not closed properly, some of the aerosol will end up in your mouth and throat. One sign of this is that the nostrils do not vibrate. In this case, it is possible that not enough saline solution or medication is reaching the paranasal sinuses.*

4) The soft palate is the extension of the hard palate.

Procedure

DANGER


Life-threatening situation due to mixing up tubes

If tubing systems that belong to other devices are close by (e.g., for infusions), there is a possible danger of mixing up different connection tubes.


- Make quite sure that the other ends of the tubes connected to the compressor are attached to the PARI nebuliser.

Proceed as follows in order to carry out the treatment:

- Take the nebuliser out of its holder and hold it upright.
- Verify that all parts are firmly connected to each other.


- Switch the compressor on.
 **DANGER! Danger of death by electrocution in the case of device fault!** Switch the compressor off and disconnect the power plug from the mains socket immediately if a fault is suspected (e.g., following a fall or if there is a smell of

burning plastic). If there is a fault in the device, it may be possible to come into contact with live parts. This in turn may lead to an electric shock.

- Check that an aerosol is being generated before you begin the treatment.
- Sit in an upright position and relax.
- Constrict one nostril.

To do this, either press one finger lightly against the side of the nose or use the nasal plug supplied.

i *Make sure that the nostril is not completely blocked (e.g. by pressing closed with a finger) because this might create a feeling of too much pressure.*

The nasal plug has a small aperture. This aperture reduces the feeling of pressure.

- Press the nasal joining piece of the nebuliser against the free nostril so that the nostril is closed tightly.
- Close your soft palate [see: Tips for closing the soft palate, page 17].

The aerosol is automatically delivered to the paranasal sinuses. Therefore, you should hold your breath during the treatment. In order to breathe, pause the treatment briefly.

After about 2.5 minutes, carry out the treatment as described for a further 2.5 minutes on the other nostril. The treatment is complete after 5 minutes in all.

i *Some residual fluid will remain in the nebuliser after the end of the treatment.*

Reducing feeling of pressure


If you feel an uncomfortable level of pressure while carrying out the therapy, reduce the pressure:

- reduce the pressure on your free nostril and/or
- enlarge the hole in the nasal plug.

To enlarge the hole, proceed as follows:

i *The larger the hole, the less the feeling of pressure. Proceed carefully to avoid making the hole too large and reducing the pressure too much. You must be able to feel the vibration all the time during the therapy.*

- Push out the inside of the nasal plug and cut off the tip at the edge indicated using a pair of scissors.
- Push the inside of the nasal plug back in.

3.3 Ending the treatment

To end the treatment, proceed as follows:

- Switch the compressor off.
- Place the nebuliser back in the holder on the compressor.
- Disconnect the power plug from the mains socket.

i *Complete disconnection from the mains is only certain when the power plug has been unplugged from the socket.*

4 REPROCESSING

The product components must be cleaned thoroughly immediately after each application and disinfected at least once a week. In the case of an acute infection, disinfection must be carried out every day.

When used by high-risk groups (e.g., patients with cystic fibrosis), the components must be disinfected at least once a day.

If the nebuliser is used **in professional environments**, follow the information on reprocessing included in the appendix at the end of these instructions for use.


4.1 Preparation

- Pull the two tubes and the adapter off the nebuliser.
- Make sure that all residual volume is removed from the nebuliser.
- Detach the adapter from the tube.
- Dismantle the nebuliser into its individual parts.

4.2 Cleaning

i *The PARI SINUS tubing system cannot be cleaned or disinfected. Regarding care of the tubing system, follow the instructions in the corresponding section [see: Care of the PARI SINUS tubing system, page 23].*

- Briefly rinse all parts used in running drinking water beforehand.
- Place all disassembled components in warm drinking water with a little dishwashing liquid for about 5 min.

- Rinse off all parts thoroughly in drinking water.
- Shake the water out of all of the parts.

4.3 Disinfecting

After cleaning, disinfect all of the **disassembled parts** (only parts that have been cleaned can be disinfected effectively). The recommended disinfection procedures are described below. Descriptions of other validated disinfection procedures are available from the manufacturer or dealer upon request.

i *The PARI SINUS tubing system cannot be cleaned or disinfected. Regarding care of the tubing system, follow the instructions in the corresponding section [see: Care of the PARI SINUS tubing system, page 23].*

CAUTION

Risk of infection due to moisture

Moisture encourages the growth of bacteria.

- Remove all parts from the pot or disinfectant as soon as the disinfection process is finished.
- Dry all parts.

In boiling water

- Place all the **individual parts** in boiling water for at least 5 minutes. Use a clean pot and drinking water.

NOTE! Risk of damage to plastic parts! Plastic will melt if it comes into contact with the hot base of the pot. Make sure there is enough water in the pot to prevent the individual parts from touching the pot base.

- Shake the water out of all of the parts.

Using a standard thermal disinfecter for baby bottles (not a microwave oven)

CAUTION

Risk of infection due to inadequate disinfection

Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection.

- Make sure that the disinfecter is clean and operating properly before every disinfection process.
- Allow the disinfection to continue until the disinfecter switches off automatically or the minimum disinfection time stated in the instructions for use of the disinfecter has elapsed. Do not switch the device off prematurely.

Use a thermal disinfecter with a runtime of at least 6 minutes. Regarding the performance of the disinfection, the duration of the disinfection procedure and the quantity of water required, follow the instructions for use of the disinfecter you are using.

4.4 Care of the PARI SINUS tubing system

Dry the PARI SINUS tubing system after every use:

- Insert the matching tubing connector on the first of the two tubes in the compressed air connector on the compressor.
- Switch the compressor on.
- Allow air to blow through the tube until the air has removed all moisture from the tube.
- Switch the compressor off.
- Disconnect the tube from the compressor, and attach the other tube to the compressed air connection.
- Repeat the process as described for the first tube.

4.5 Inspecting

Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discoloured parts.

4.6 Drying

After each cleaning and disinfection, place all product components on a dry, clean and absorbent surface and let them dry completely.

4.7 Storage

Store this product as described below:

- Wrap all individual components in a clean, lint-free cloth (e.g., a tea towel).
- Store this product in a dry, dust-free place.

5 TROUBLESHOOTING

Fault	Possible cause	Remedy
No aerosol is coming out of the nebuliser.	The nebuliser nozzle insert is blocked.	Clean the nebuliser.
	The compressed air supply tube is not connected correctly to the compressed air connection on the compressor or the nebuliser lower section.	Check that the tubing connectors are connected firmly to the compressor and the nebuliser.
	The tube attached to the compressed air connection on the compressor is leaking.	Replace the defective tube or the entire PARI SINUS tubing system.

6 TECHNICAL DATA

6.1 General nebuliser data

Size ⁵	10 cm × 10 cm × 4 cm
Weight ⁵	31 g to 33 g
Operating gases	Air
Minimum compressor flow	3.0 l/min.
Minimum operating pressure	0.5 bar / 50 kPa
Maximum compressor flow	6.0 l/min.
Maximum operating pressure	2.0 bar / 200 kPa
Minimum fill volume	2 ml
Maximum fill volume	8 ml

6.2 Aerosol data

The aerosol characteristics of the PARI LC SPRINT SINUS nebuliser were determined using the PARI SINUS2 compressor. The particle size was measured with the “Next Generation Impactor” (NGI). The measurement was carried out with an ambient temperature of 23 °C and 50% relative humidity. 2 ml salbutamol (1 mg/1 ml) were nebulised. The Jet Flow was measured to be 4.6 l/min.

Nozzle insert (orange)	Nominal compressor flow (4.6 l/min – 1.5 bar)
MMAD [μm] ⁶	3.0
GSD ⁷	2.18
Aerosol fraction [% > 5 μm]	66.6
Aerosol fraction [% < 2 μm]	28.7

5) Without bend, nasal joining piece and adapter; unfilled.

6) MMAD = Mass Median Aerodynamic Diameter

7) GSD = Geometric Standard Deviation

Nozzle insert (orange)	Nominal compressor flow (4.6 l/min – 1.5 bar)
Aerosol fraction [% > 2 µm < 5 µm]	38.0
Aerosol fraction [% > 5 µm]	33.4
Aerosol output [ml]	0.78
Aerosol output rate [ml/min]	0.15
Residual volume [ml] (gravimetric)	0.95
Percentage of fill volume emitted per minute [%/min]	7.54

7 MISCELLANEOUS

7.1 Disposal

All product components can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

7.2 Labelling

The following symbols can be found on the product and/or the packaging:

	Legal manufacturer
	The product satisfies the basic requirements as set forth in Appendix I of Directive 93/42/EEC concerning medical devices.
	Consult instructions for use
	Item number
	Production batch number, lot

APPENDIX: Reprocessing in professional environments (for only one patient)

Nebuliser

The following overview of the processing steps in professional environments applies to the following products:

- Nebuliser (all components except the PARI SINUS tubing system)

1. Preparation

Disassemble the product [see: Preparation, page 21].

Check:

- Expiry date of the cleaning agent/disinfectant
- Processing limits reached?

2. Cleaning and disinfection

Manual cleaning:	pH-neutral, enzymatic cleaner, e.g., Korsolex® Endo Cleaner (Bode) or Bodedex® forte (Bode) Use: In accordance with manufacturer information, if dirt is visible use a brush if necessary.
Manual disinfection:	<p>With aldehyde-containing instrument disinfectant, e.g., Korsolex® Basic (Bode) Use: In accordance with manufacturer information Active agent basis: – Aldehyde donor ¹ – Aldehyde ¹</p> <p>With aldehyde-free instrument disinfectant Use: In accordance with manufacturer information Active agent basis: Quaternary ammonium compound <i>Info: Bomix® plus cleans and disinfects in a single work step</i></p>
Mechanical cleaning with disinfection:	Neutral cleaning agent, e.g., neodisher® Medizym (Dr. Weigert) or alkaline cleaning agent, e.g., neodisher® MediClean forte 0.5% (Dr. Weigert) in conjunction with neutraliser, e.g., neodisher® Z (Dr. Weigert) Equipment: Cleaning device and disinfectant in conformance with DIN EN ISO 15883, e.g., RDG G7836 CD (Miele) Vario TD programme or comparable valid programmes

1) Additional activity spectra: tuberculocidal, mycobactericidal, fungicidal
According to the information supplied by the disinfectant manufacturers, the agents listed are effective against the gram-positive bacteria *Staphylococcus aureus* and *Enterococcus hirae* and against the gram-negative bacteria *Escherichia coli*, *Pseudomonas aeruginosa* and *Proteus mirabilis* as well as the yeast-like fungus *Candida albicans*.

3. Steam sterilisation

Equipment:

- Steam steriliser (preferably with fractionated pre-vacuum) in accordance with DIN EN 285 or DIN EN 13060 (Type B)
- Sterile barrier system in accordance with DIN EN 11607

Temperature / Duration:

134 °C for at least 3 min.

4. Visual inspection & storage

Check:

Inspect all individual parts. Replace any broken, misshapen or seriously discoloured parts.

Storage location:

- dry
- dust-free
- protected from sources of contamination
optional: Use sterile packaging

Processing limits

max. 300 disinfection cycles, max. 1 year

PARI SINUS tubing system

1. Preparation

Check the product:

- Expiry date of the cleaning agent/disinfectant
- Processing limits reached?

2. Cleaning and disinfection

Manual cleaning:	not applicable
Manual disinfection:	not applicable

<p>Mechanical cleaning with disinfection:</p>	<p>Alkaline cleaning agent, e.g., neodisher® MediClean forte (Dr. Weigert) in conjunction with neutralising agent, e.g., neodisher® Z (Dr. Weigert)</p> <p>Equipment:</p> <ul style="list-style-type: none">– Cleaning device and disinfector in conformance with DIN EN ISO 15883, e.g., RDG G7836 CD (Miele)– Special baskets for Miele instrument dishwasher– Compressed air source for blowing dry <p>Vario TD programme or comparable valid programmes</p>
---	---

3. Steam sterilisation

not usable

4. Visual inspection & storage

Inspect all individual parts. Replace any broken, misshapen or seriously discoloured parts.

Storage location:

- dry
 - dust-free
 - protected from sources of contamination
- optional: Use sterile packaging

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