

Wofasept® FL spray disinfection

Ready-to-use surface disinfectant



For odourless disinfection of washable surfaces in the working area of medical practices, hospitals and other medical facilities as well as laboratories, food industry and gastronomy. Can be used when alcoholic disinfection is undesirable.

New:
Bottle with
practical foam
and spray
function



Highly effective microbicide

Wofasept® FL spray disinfection effectively kills bacteria, funguses, lipophilic viruses (including HBV/HIV) and human rotaviruses



No protein related faults

Insensitive to contamination such as blood, protein, stool, sputum, pus, etc.



No allergic potential

Wofasept® FL spray disinfection is free of chlorine, aldehyde and phenol and does not contain any substances from which an allergic potential is to be expected.



Application and dosage:

Spray or foam undiluted wipe and wipe away with a dry tissue after application, if required.

	Application time
for surface disinfection in hospitals, practices and other medical facilities	
• kills bacteria and funguses	5 min
• kills hepatitis B and HIV	1 min
• kills herpes viruses	15 min
• kills human rotaviruses	1 min
• kills vacciniaviruses	5 min
• kills noroviruses	30 min
• kills influenza A-viruses	1 min
Disinfects solariums and washable surfaces	5 min
for food industry	30 min



Wofasept® FL spray disinfection

Ready-to-use surface disinfectant



Composition:

0,24 % Alkyl(C12-C16)dimethylbenzylammoniumchlorid, 0,24 % Alkyl(C12-C14)dimethyl(ethylbenzyl)ammoniumchlorid, 0,24 % Didecyldimethylammoniumchlorid

Possible interactions:

The combined use of Wofasept® FL spray disinfection with other disinfectants or cleaning agents should not be used, as chemical interaction can lead to a considerable loss of effectiveness. Thoroughly rinsing with water before disinfection is necessary important after using anionic cleaning agents.

Material compatibility:

When using it on sensitive surfaces and materials, a preliminary test with the intended concentration and exposure time at an inconspicuous place is recommended.

Hazard statements:

H319 Causes serious eye irritation .

Precautionary statements:

P264 Wash skin thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P337 + P313 If eye irritation persists get medical advice/attention. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. Remove any existing contact lenses if possible. Continue rinsing.

Rinse children's toys with water after disinfection.

Storage:

Keep out of the reach of children. Store at room temperature in closed original package, protected from direct sunlight, heat and moisture in a well-ventilated place.

Disposal:

Add small amounts diluted with plenty of water to the wastewater.

Larger quantities must be disposed of in accordance with local official regulations via an approved waste disposal company.

Stability:

If stored properly, it can be kept for at least 5 years without any loss of effectiveness.

Use biocide safely. Always read labelling and product information before use.

Biocide registration:

BAuA-Reg-Nr.:
[German Federal Institute for Occupational Safety and Health]

N-22126 (product type 2)
N-76931 (product type 3)
N-22127 (product type 4)

package sizes:

500 ml	spray bottle	PZN: 06878235
5 l	can	PZN: 06878229



Your partner for professional hygiene management.

KESLA HYGIENE AG

Keslastraße 2
06803 Bitterfeld-Wolfen

Tel.: +49 (0) 3494 - 6995 0
Fax: +49 (0) 3494 - 6995 29
www.kesla.de
info@kesla.de



DIN EN ISO 9001

No guarantee can be given for the general validity of the information given in this product information. All information reflect our current level of knowledge and experience, and do not warrant any legal liability.

Updating: PI 08.20